Cotilos tipo Müller cementados
Cemented cups Müller type

surgival
GRUPO COSÍAS

Características técnicas Technical details

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Surgival Müller type cups are characterised for being hemispherical cups manufactured entirely from ultra-high-molecular-weight polyethylene (UHMWPE). The hemispherical design enables the maximum bone capital to be conserved, guaranteeing uniform load distribution to the bone thanks to a narrow contact with the acetabular bone of less than 0.5 mm.

They are cemented cups, which means that on their convex external surface they have a special grooved design aimed at improving the cup's primary fastening. This primary fastening takes place on the cup's impact on the acetabular cavity, which is previously filled with bone cement.

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This impaction ensures the bone cement permeates the mechanised grooves in the cup and that, once the bone cement has set, they act as retainers to prevent the cup from rotating or loosening. The purpose of these grooves is also to increase the specific surface on contact with the acetabulum, providing greater peripheral contact and limiting the appearance of possible micro-movements that may cause prosthetic loosening.

On the outside, all cups in the range have the same finishes in terms of grooves and channelling, depending on the number and placement of the cup's acetabular size or diameter.

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The range is made up of the variants Antiluxation Cup and Standard.

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All the cups have a marker ring whose purpose is to assist the surgeon in checking the correct acetabular placement of the implanted cup. Since UHMWPE is a radiotransparent material, this opaque-radio marker ring is necessary for later post-operative clinical monitoring.

The three external orifices are for lodging the specific instruments to facilitate positioning, orientation and impactation on the acetabulum.

The minimum polyethylene thickness (UHMWPE) in the range of Surgival Müller type cups is greater than the minimum value specified in the norm. This greater thickness provides greater dispersion of pressures between the insert/cement/bone interfaces, resulting in greater durability of the cup and the fastening.

Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5

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1. Acetabulum exposure and drilling

To correctly prepare the acetabulum (elimination of all the soft parts, extirpation of the articular cartilage and making perforations in the acetabular bone), there must be a complete and wide view of this cavity. Placing Hohmann separators on the anterior and posterior edges and below the acetabular groove obtains a complete view of this cavity.

Use progressive sized milling cutters to eliminate the acetabular cartilage until there is a bleeding subchondral bone. (Fig.1)

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2. Test cups

Assemble the tester in its sleeve and place it in the acetabulum to check its adaptation. The design of the testers enables the base of the cup cavity to be viewed to confirm it is adapted to the acetabular surface. (Fig. 2)

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3. Acetabulum perforations

With the help of a drill bit with a safety stop, make several perforations on the acetabular bone surface. Next insert the cement in these perforations. (Fig. 3)

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4. Placing the acetabular cup

After meticulously washing and cleaning the acetabulum, insert a portion of the acrylic cement in a paste state at pressure into the cavity. The polyethylene acetabular cup is mounted in the guide-impactor and inserted with the most suitable inclination and anteversion, pressing against the bone wall of the acetabulum to express the excess cement. Keep the guide-impactor with the acetabular cup immobile until the cement has set. (Fig. 4)

The implant selected must have a diameter of 2 or 4 mm less than the acetabular cavity in order to obtain a sufficient thickness of cement layer around it.

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NOTE:

Pressing the button located on the impactor handle assists with its separation from the implanted cup thanks to the mechanism, which helps with expulsion while keeping the cup fixed in the acetabular cavity. (fig 5)

Instrumental Instruments

Set completo instrumental cotilo antiluxante
Antiluxation cup instrument complete set

Ref. A1512000

Bandeja instrumental cotilo antiluxante
Antiluxation cup instrument tray

Ref. A1512001

Caja vacía 1 filtro 60x30x11 cm
Empty container 1 filter 60x30x11 cm

Ref. B4500001

1	Fresa acetabular Acetabular reamer
Ref. A1701044	44 mm
Ref. A1701046	46 mm
Ref. A1701048	48 mm
Ref. A1701050	50 mm
Ref. A1701052	52 mm
Ref. A1701054	54 mm
Ref. A1701056	56 mm
Ref. A1701058	58 mm

2	Cotilos de prueba Trial cups
Ref. A2500844	44 mm
Ref. A2500846	46 mm
Ref. A2500848	48 mm
Ref. A2500850	50 mm
Ref. A2500852	52 mm
Ref. A2500854	54 mm
Ref. A2500856	56 mm
Ref. A2500858	58 mm

3	Mango porta-fresa acetabular Handles for acetabular reamers
Ref. A1701075	2 unid. / 2 units

4	Fresa para perforaciones acetabulares Reamers for acetabular drillings
Ref. F1711000	

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6	Mango impactor cotilo antiluxante Antiluxation cup impactor sleeve
Ref. A1512500	

7	Acople universal Universal coupling
Ref. A1700380	

Implantes Implants

Cotilos tipo Müller
Müller type cups

Antiluxante / Antiluxation

Ref. A1512244E	28 x 44
Ref. A1512246E	28 x 46
Ref. A1512248E	28 x 48
Ref. A1512250E	28 x 50
Ref. A1512252E	28 x 52
Ref. A1512254E	28 x 54
Ref. A1512256E	28 x 56
Ref. A1512258E	28 x 58

Standard *

Ref. A1511244E	28 x 44
Ref. A1511246E	28 x 46
Ref. A1511248E	28 x 48
Ref. A1511250E	28 x 50
Ref. A1511252E	28 x 52
Ref. A1511254E	28 x 54
Ref. A1511256E	28 x 56
Ref. A1511258E	28 x 58

Nota: Existe un mango específico para la colocación del Standard. Consulte a su distribuidor o fabricante.
Note: There is a specific handle for the Standard. Consult your distributor or manufacturer.

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