


OSTEOMED

A TRADITION OF QUALITY AND INNOVATION

OSTEOHARVESTER™

Harvester • Morselizer • Collector • Delivery System



MORSELIZED AUTOGENOUS BONE... IN A SINGLE AUTOMATED STEP.

- CAN ALSO BE USED AS AN IN-LINE SUCTION BONE COLLECTOR
- MINIMALLY INVASIVE METHOD FOR HARVESTING CORTICAL AND CANCELLOUS BONE PARTICULATE
- COMPATIBLE WITH ANY CONTRA-ANGLE HANDPIECE
- AFFORDABLE SYSTEM WITH MINIMAL COST DISPOSABLES



PRODUCTS THAT IMPROVE PEOPLE'S LIVES

OSTEOHARVESTER™

Harvester • Morselizer • Collector • Delivery System

BONE HARVESTER (ITEMS NOT ACTUAL SIZE)

REUSABLE COMPONENTS




DISPOSABLE COMPONENTS



BONE COLLECTOR (ITEMS NOT ACTUAL SIZE)

REUSABLE COMPONENTS



DISPOSABLE COMPONENTS



Can be used with:
OSTEOPower™
SURGICAL HANDPIECE SYSTEM

PART NUMBER	DESCRIPTION
230-0001	RESERVOIR
230-0002	HARVEST CAP
230-0003	COLLECTOR CAP
230-0005	PLUNGER
230-0017	PLUG
230-0018	SUCTION ADAPTER, 1/2" DIAMETER
230-0019	SUCTION ADAPTER, 1/4" DIAMETER
230-0021	COLLECTOR WAND
230-0029	DISH
230-0032	CONTRA-ANGLE DISPOSABLE SET, 7MM PLUNGE, 3 PC
230-0034	CONTRA-ANGLE DISPOSABLE SET, 11MM PLUNGE, 3 PC
230-0036	COLLECTOR FILTER



3885 Arapaho Road
Addison, Texas 75001
Customer Service: 800.456.7779
Main: 972.677.4600
Fax: 800.390.2620
www.osteomedcorp.com

P/N 030-1193 Rev.G